Borescope Imaging for Surgical Safe Path Planning

Xingchi He, Saumya Gurbani, Alperen Degirmenci

Mentors: Dr. Taylor, Dr. Iordachita

Clinical Advisor: Dr. Chien








Outline

- Background & Motivation
- Project Goal
- Technical Approach
- Deliverables
- Timeline
- Dependencies
- Reserposibilities
- Management Plan
- Reading List





Background & Motivation








Project Goal

- The goal of this project is to create safe insertion paths for surgeons using a borescope.
- Components of this project:
 - A hardware adapter which mounts the borescope system on the EyeRobot
 - Software to generate safe paths
 - Software to register the safe path / virtual fixtures with subsequenct surgical tools.





Technical Approach / Work Flow







Deliverables

	Deliverable	Status				
Minimum	Hardware to mount borescope system on to the EyeRobot setup.	EXPECTED 04/29/2011				
	Develop software for "manual" safe-path generation	EXPECTED 05/03/2011				
Expected	Develop software to register the boroscope & safe path with subsequent surgical tools	EXPECTED 05/17/2011				
Maximum	Generate virtual fixtures from safe path	EXPECTED SUMMER 2011				
	Automate the safe path generation process	EXPECTED SUMMER 2011				
	Use optical flow to create a 3D reconstructions of the cochlea	EXPECTED SUMMER 2011				







Borescope Safe Path Planning Timeline (4/5/11) FERC | CISST

STAGE	WEEK TASK	8	9	10	11	12	13	14	Summer				
0	Project Planning												
1	Build Light Source/Cable adapters												
2A	First prototype of mounting adapter												
3	Software for "manual" safe path creation												
2B	2 nd prototype of mounting adapter												
4	Borescope-Surgical Tool Registration												
	Final Presentation + Demo												
5	Automated Safe Path Creation												
6	3D Reconstruction Using Optical Flow												
7	Virtual Fixture Generation												
© 2011 LCSR – ERC CISST Johns Hopkins University		April		May		July/August							

Dependencies

Dependency	Status as of 04/05/2011					
Micro-Borescope	DAMAGED (resolution: use backlighting)					
Camera for recording boroscope images	RESOLVED 04/05/2011					
Light Source/Cable (Purchased or Created)	RESOLVED 04/08/2011					
Access to EyeRobot and Marcin	RESOLVED 02/28/2011					
Access to Machine Shop and Rapid Prototyping Machine	RESOLVED 02/14/2011					





Management Plan

- 1. Weekly meeting with Dr. Iordachita(if he is available)
- 2. Weekly meetings with each other on Tuesday and Thursday
- 3. 50 man hrs total per week





Reading List

- [1] Cohen, Noel L. and J. Thomas Roland Jr. "Complications of Cochlear Implant Surgery". Cochlear Implans. Ed. Susan B. Waltzman and J. Thomas Roland. Thieme Medical Pub, 2006. 126-132.
- [2] I. Fleming, M. Balicki, J. Koo, I. Iordachita, B. Mitchell, J. Handa, G. Hager, and R. Taylor, "Cooperative robot assistant for retinal microsurgery.," *Medical image computing and computer-assisted intervention : MICCAI ... International Conference on Medical Image Computing and Computer-Assisted Intervention*, vol. 11, Jan. 2008, pp. 543-50.
- [3] C.G. Wright, P.S. Roland, and J. Kuzma, "Advanced bionics thin lateral and Helix II electrodes: a temporal bone study.," *The Laryngoscope*, vol. 115, Nov. 2005, pp. 2041-5.
- [4] D. Zhang et. al., "Inspecting the cochlear scala tympanic with flexible and semi-flexible micro-endoscope." Journal of Clinical Otorhinolaryngology Head and Neck Surgery, vol. 20 issue 4, 2006. pp 169-71.



