Intracardiac Mapping of Clinical Markers for Prediction of Ventricular Tachyarrhythmia

Computer Integrated Surgery II
Spring Semester, 2012
Markus Kowalsky and Sindhoora Murthy

Introduction

- Ventricular tachyarrhythmias play a key role in sudden cardiac death which claims the lives of approximately 400,000 people per year in the United States.
- One of the current therapies for preventing tachyarrhythmias, an implantable cardioverter-defibrillator, is only effective in ~65% of patients
- Recently, Dr. Larisa Tereshchenko, showed that the sum absolute QRST integral is an accurate predictor for ventricular tachyarrhythmia
- Unfortunately, the current method for mapping the electrical activity in the heart involves a catheter-fed electrode and is unable to map individual heart beats
- Our project seeks to allow physicians to understand the electrical activity in the heart through noninvasive 120-lead ECGs

Methods


- ECG Fiducial Points were determined semi-automatically using a modified version of the algorithm detailed in "A QT Interval Detection Algorithm Based on ECG Curve Length Transform" (Zong et al. 2006)
- From these detected fiducial points, various clinical markers were mapped to their corresponding nodes in a 120/291 node heart-torso model. The remaining surface was calculated through interpolation.

Example Method

Curve Length Data with Q and S Markup Original Data with Fiducial Points 2500 - ECG * Q 2000 * R * s Length Transform Crossing Point (tc) tc-125ms min tc+125ms max 1000 Sample Number Zeroed QRS Data Zeroed QRS Curve Length Data with End of T Markup 1200 800 Curve Length Data Secant Line End of Line Start of Line Zoomed In Zeroed QRS Curve Length Data with End of T Markup Rotated T-end detection 1595 1590 1592 Curve Length Data End of Line Start of Line


Heart-Torso Model

Rotated Curve Length Data




ا 1570 ك

Ο 1565



Results

- Areas with elevated SAI QRST and/or high peak-to-peak voltage matched closely with hypertrophied tissue as seen in MRI scan
- QT interval map shows dispersion of depolarization at the edge of arrhythmogenic scar tissue
- SAI QRST may be a surrogate for peak-to-peak voltage
- The placement of each orthogonal lead for traditional SAI QRST measurement is very important



Future Work

- Additional Validation with 6 additional patients
- Explore whether adding additional information (tissue conductivity, edema status, etc.) to model improves inverse solution and intracardiac mapping
- Further study of these maps to determine optimal placement of left ventricular lead in bi-ventricular pacing

Acknowledgements

We would like to thank Dr. Larisa Tereshchenko for her support throughout this semester's project and Dr. Fady Dawoud for providing all of the necessary data and guidance with construction of the maps.



