


NSF Engineering Research Center  
for Computer Integrated Surgical  
Systems and Technology



# MICRON RANGE-OF-MOTION VISUALIZATION

Check Point Presentation II

Team-14

## **Preetham Chalasani**

**Department of Computer Science**

**The Johns Hopkins University**

**pchalas1@jhu.edu**

**Mentors - Dr.Russell Taylor, Marcin Balicki, Balazs Vagvolgyi**

**WHITING  
SCHOOL OF  
ENGINEERING**

THE JOHNS HOPKINS UNIVERSITY







# CONTENTS



# SUMMARY





- Need :- Surgeons don't always know the position of the micron in its range of motion
- Goal :- Develop a visual alert assistance system for the surgeons dealing with very small anatomy.




# MICRON

- Tremor Cancellation
- Move actively to compensate



# SOLUTION






# TIMELINE

	Week Starting with													
	Feb. 4	Feb. 11	Feb. 18	Feb. 25	Mar. 4	Mar. 11	Mar. 18	Mar. 25	Apr. 1	Apr. 8	Apr. 15	Apr. 22	Apr. 29	May 6
PHASE - I														
Understanding CISST and StereoVision libraries														
Setting up development Environment														
Understanding the Existing Framework														
Create a test Application														
Include some overlays														
Develop Application using simulated data														
Communicate with the micron and get the information														
Develop Application using Micron data														
Integrate Tool Tracker														
Continuous Feedback														
Rigorous Testing														
Debugging														
Include the micron tip workspace														
Improve the tracker														





# DELIVERABLES



# CURRENT PROGRESS

Component Connections

Random Data

Actual Data


Experimental Data

Audio Feedback

Micron Range (Cube) Display



# INTERNAL COMPONENT CONNECTION



Source : <https://trac.lcsr.jhu.edu/cisst/wiki/cisstMultiTaskTutorial>

# MY CONNECTIONS

## Provided

devMicron

micronLimitsBehavior

devMicron

devOpenAL



## Required

micronLimitsBehavior


micronPainter

AsapGUI


micronLimitsBehavior



# GUI FLOW DIAGRAM



# FLOW DIAGRAM




# EXPERIMENTAL DATA

- Calibration (ASAP)
  - Yaw : 22.5
  - Pitch : 50
  - Roll :-10
- Range – Cylindrical
  - Height : 4mm
  - Radius : 2mm
- 1mm : 100 px



# EXPERIMENTAL FEATURES



Summary

Background

Progress

Deliverables

Dependencies

Timeline

Reading List



# AUDIO FEEDBACK


- **Safe**
  - $(\text{Tip} - \text{HomeTip}) < (\text{MicronRadius} - 15 \text{ px})$
- **Warning**
  - $(\text{MicronRadius} - 15 \text{ px}) < (\text{Tip} - \text{HomeTip}) < (\text{MicronRadius} - 8 \text{ px})$
- **Red Alert**
  - $(\text{Tip} - \text{HomeTip}) > (\text{MicronRadius} - 8 \text{ px})$





# AUDIO FEEDBACK

Safe



Micron Tip


Micron Workspace





# AUDIO FEEDBACK

Warning



● Micron Tip


— Micron Workspace





# AUDIO FEEDBACK

Alert





● Micron Tip

— Micron Workspace



# MICRON RANGE DISPLAY



Without Depth

Summary

Background

Progress



Deliverables

Dependencies

Timeline

Reading List

# MICRON RANGE DISPLAY



Without Depth

Summary

Background

Progress



Deliverables

Dependencies

Timeline



Reading List

# MICRON RANGE DISPLAY



## With Depth

# VIDEO



# CONFUSION !!

Cube

?????

Sphere




# PROBLEMS FACED

- Access to Micron
- Build/Compile Delay (Lot of Sleeps)
- 3D on 2D



# 3D ON 2D



Summary

Background

Progress

Deliverables


Dependencies

Timeline

Reading List



# UPDATED DELIVERABLES





# DEPENDENCIES

Dependency	Source	Status/Comments	What If ??
PC or Laptop	Self	Acquired	Project Delayed
CISST and Stereo Vision Libraries	Open Source-Online	Installed	Custom Libraries
QT Creator - IDE	Open Source-Online	Installed	Use other free IDEs available
Material to understand Micron better	Dr.Russel Taylor	Acquired	Learn Myself
Documentation of previous work	Marcin Balicki/Balazs Vagvolgyi	Acquired	Learn myself
Access to micron	Marcin Balicki/Balazs Vagvolgyi	Acquired - Not always available	Work on simulated data/Project Delayed
Access to Stereo video Microscope	Marcin Balicki/Balazs Vagvolgyi	Acquired - Not always available	Work on simulated data/Project Delayed





# UPDATED TIMELINE

	Week Starting with													
	Feb. 4	Feb. 11	Feb. 18	Feb. 25	Mar. 4	Mar. 11	Mar. 18	Mar. 25	Apr. 1	Apr. 8	Apr. 15	Apr. 22	Apr. 29	May 6
PHASE - I														
Understanding CISST and StereoVision libraries														
Setting up development Environment														
Understanding the Existing Framework														
Create a test Application														
Include some overlays														
Develop Application using simulated data														
Communicate with the micron and get the information														
Develop Application using Micron data														
Integrate Tool Tracker														
Continuous Feedback														
Rigorous Testing														
Debugging														
Include the micron tip workspace														
Improve the tracker														





# READING LISTS

- [1] B. C. Becker, S. Voros, R. A. MacLachlan, G. D. Hager, and C. N. Riviere, “Active Guidance of a Handheld Micromanipulator using Visual Servoing”, in IEEE International Conference on Robotics and Automation, Kobe, Japan, May 12-17, 2009. pp. 339-344.
- [2] B. Becker, R. MacLachlan, and C. Riviere, “State estimation and feedforward tremor suppression for a handheld micromanipulator with a Kalman filter”, in EEE RSJ Int Conf Intell Robot Syst, 2011. pp. 5160-5165. NIHMSID: 345014.
- [3] B. Becker, R. MacLachlan, L. Lobes, and C. Riviere, “Vision-Based Retinal Membrane Peeling with a Handheld Robot”, in IEEE Int Conf Robot Autom, 2012. pp. 1075-1080. NIHMSID: 368417.
- [4] B. Becker, S. Yang, R. MacLachlan, and C. Riviere, “Towards vision-based control of a handheld micromanipulator for retinal cannulation in an eyeball phantom”, in Proc IEEE RAS EMBS Int Conf Biomed Robot Biomechatron, 2012. p. accepted for publication. NIHMSID: 368431.
- [5] B. Gonenc, M. A. Balicki, J. Handa, P. Gehlbach, C. N. Riviere, R. H. Taylor, and I. Iordachita, “Preliminary Evaluation of a Micro-Force Sensing Handheld Robot for Vitreoretinal Surgery”, in IEEE/RSJ International Conference on Intelligent Robots and Systems (IROS), Vilamoura, Algarve, Portugal, 7-12 October, 2012. pp. 4125-4130.





# READING LISTS

- [6] R. MacLachlan, B. Becker, J. Cuevas-Tabarés, G. Podnar, L. Lobes, and C. Riviere, "Micron: an actively stabilized handheld tool for microsurgery", IEEE Trans Robot., vol. 28- 1, pp. 195-212, 2012. NIHMSID:345015.
- [7] S. Yang, M. Balicki, R. A. MacLachlan, X. Liu, J. U. Kang, R. H. Taylor, and C. N. Riviere, "Optical Coherence Tomography Scanning with a Handheld Vitreoretinal Micromanipulator ", in IEEE Engineering in Medicine and Biology Conf, San Diego, Aug 28-Sep 1, 2012. pp. 948-951. NIHMSID: 383510.
- [8] S. Yang, R. MacLachlan, and C. Riviere, "Design and analysis of 6 DOF handheld micromanipulator", in Proc IEEE Int Conf Robot Autom., St. Paul, MN, May 14-18, 2012. pp. 1946-51. NIHMSID: 368427.
- [9] B. Becker, R. MacLachlan, L. Lobes, G. Hager, and C. Riviere, "Vision-Based Control of a Handheld Surgical Micromanipulator with Virtual Fixtures", IEEE Transactions on Robotics, pp. Accepted Nov 27, 2012, 2013. NIHMSID: 429749.
- [10] M. Balicki, J.-H. Han, I. Iordachita, P. Gehlbach, J. Handa, R. H. Taylor, and J. Kang, "Single Fiber Optical Coherence Tomography Microsurgical Instruments for Computer and Robot-Assisted Retinal Surgery", in Medical Image Computing and Computer Assisted Surgery (MICCAI 2009), London, September 20-24, 2009. pp. 108-115. PMID: 20425977



# QUESTIONS?